


# Unravel the Chloroplast Genome



## *Hereditas* *Photosynthetica*

- Harvest 10 g of spinach seedlings (kept in dark for 48 hrs to reduce starch levels). Wash under tap water. Cut into 1 cm strips (directly into ice-cold buffer).
- Grind in ice-cold mortar with 20 ml of cold complete buffer for 2-4 minutes. Don't grind excessively: Homogenate should be "chunky" (with some incompletely homogenized material).
- Filter homogenate through 50  $\mu\text{m}$  nylon mesh into a 40 ml Oakridge centrifuge tube. A filter funnel will avoid messes. The material in the nylon mesh can be gently squeezed to release more chloroplasts into the tube.
- Centrifuge for 20 minutes at 1 000 rpm at 4°C. This will pellet nuclei and cell debris, leaving mostly chloroplasts (intact and broken) in the supernatant.
- Gently pour supernatant into a new Oakridge centrifuge tube and spin at 4 000 rpm for 20 minutes at 4°C. This will pellet the chloroplasts.
- Discard supernatant and resuspend pellet (gently!) in 200  $\mu\text{l}$  ST buffer.
- Transfer sample to microfuge tube and add equal volume of equilibrated phenol. Shake intensively for 5 minutes, centrifuge 5 minutes at 12 000 rpm. There will be 2 phases, green on bottom and clear on top.
- Transfer the clear upper layer to a new microfuge tube. Add equal volume of 1:1 phenol:chloroform. Shake 5 minutes and centrifuge 5 minutes at 12 000 rpm.
- Transfer upper layer to new microfuge tube and add equal volume of 24:1 chloroform:isoamyl alcohol. Shake 5 minutes and centrifuge 5 minutes at 12 000 rpm.
- Transfer upper layer to new microfuge tube and add 1/10th the volume of 5 M ammonium acetate and an equal volume of isopropanol. Store at -20°C for at least 2 hours. (Until next lab session).
- Centrifuge 5 minutes at 12 000 rpm. Discard supernatant. You may see a small white pellet (this is the DNA).
- Resuspend pellet in 200  $\mu\text{l}$  of 70% ethanol.
- Centrifuge 5 minutes at 12 000 rpm.
- Discard supernatant and air-dry the pellet, ensuring no ethanol remains in the tube.
- Resuspend pellet in 10  $\mu\text{l}$  TE buffer. This is your DNA sample!
- Mix 5  $\mu\text{l}$  of DNA sample with 1  $\mu\text{l}$  of 6X loading buffer.
- Run sample on 1% agarose gel at 100 Volts for ~45 minutes.

## Photosynthesis SC/BIOL 4160